

# UltraSlim UV Transilluminator

As small as it can be

- **Smaller:**  
Single-handedly take
- **Accurate:**  
<5% CV across the viewing area
- **Sensitivity:**  
Up to 0.5 ng nucleic acid Higher S/N ratio
- **Easy:**  
Easy for observation and gel cutting



All of the UV Transilluminators are equipped with a new-design UV blocking cover which is easily detachable for access to the work surface.



Order Information and Specification:

Products	UltraSlim UV
Cat No.	UUV-01
Outer dimensions (mm)	247D x 132W x 48.5H
View size (mm)	150D x 80W
Wavelength (nm)	302 nm
Light source	UV lamp X1
Uniformity	+/- 5 %
Power	DC12V/1.5A
Weight (Kg)	1.4 kg

