

MiniCute Imager

GEL DOCUMENTATION SYSTEM

MiniCute imager, compact size, compatible with MaestroGen's UV and blue LED transilluminators.

- 100% fit to all MaestroGen's LED or UV illuminators
- Branded digital camera, 100% dark room for photography

MiniCute imager is an easy gel image system, recommend to use with UltraSlim® LED Illuminator, UltraBright LED Transilluminator and UltraBright UV Transilluminator



COMPACT

SMALL

LIGHT

EASE

Specification


Cat no.: MC-02
 Dimensions (mm): 275 (W)* 238 (D)* 250 (H)
 Filter: Amber filter and UV-cut filter
 Weight (Kg): 1.1

Camera: Canon M-series type
 Effectivepixel: 24.2 Megapixels
 LCDSize: 3" TFT Color
 Lens: 15mm ~45mm lens, F1.8 (W) – F11.0 (T)
 Mode: (P),(A),(S),(M)
 ISO: Auto,100,200~51200
 (recommend under ISO800 for operation)

Shutter Speed(s): 1/4000~30second
 Video Output: NTSC/PAL selectable (HDMI)
 Interface: USB(v.2.0),Wi-Fi
 Storage Media: SD,SDHCmemorycard
 Storage Format: JPEG(ExifVer2.3),RAW
 Dimensions(mm): 108* 66.6* 35




Canon EOS M200



UltraBright LED Transilluminator (LB-16, LB-18)



UltraSlim® Blue LED illuminator (SLB-01W)



UltraBright UV Transilluminator (MLB-21, MLB-16, MB-21, MB-16)