


DNA Ladder

Accurate size and mass estimations

DNA Ladders comprising of fragments of specific base pair length are available for estimation of molecular weights from 100bp to 10Kb. For easy determination of band size, several DNA Ladders include bands that are two to three times more intense.



Order Information and Specification:

Products	AccuRuler 100 bp DNA Ladder	AccuRuler 100 bp Plus DNA Ladder	AccuRuler 1Kb Plus DNA Ladder
Cat No.	02001-500	02002-500	02004-500
Size	50µg / 500µl	50µg / 500µl	50µg / 500µl
Range	100 bp to 1500 bp. The 500 bp and 1,500 bp bands have increased intensity	100 bp to 3,000 bp. The 500 bp and 1,500 bp bands have increased intensity	100 bp to 10 kb. The 1 kb and 3kb bands have increased intensity
Mass/band	0.5 µg/ 5 µl	0.5 µg/ 5 µl	0.5 µg/ 5 µl