


Protein Ladder

Sharp and stable

Protein Ladders comprising of fragments of specific base pair length are available for estimation of molecular weights from 10kDa to 180kDa.



Order Information and Specification:

Products	AccuRuler RGB Prestained Protein Ladder	AccuRuler RGB Plus Prestained Protein Ladder
Cat No.	02101-250	02102-250
Size	250µl, reference bands – red at ~75 kDa and green at ~25 kDa.	250µl, reference bands – red at ~75 kDa and green at ~25 kDa.
Range	10 to 180 kDa	8 to 245 kDa
Mass/band	0.5 µg/ 5 µl	0.5 µg/ 5 µl