

Product name: AccuRuler RGB Prestained Protein Ladder

Data updated: 01/13/2022

Catalog number : 02101-250

 MaestroGen Inc.
 3F, No.18, Ln. 568, Niupu E. Rd.,
 Xiangshan Dist., Hsinchu 30091, Taiwan
 +886-3-5308665

1. IDENTIFICATION OF THE SUBSTANCE

Application of the substance For laboratory use
 the preparation: Suitable for use as molecular weight standards for SDS-PAGE electrophoresis

Manufacturer/Supplier: MaestroGen ,Inc.,
 Address: 3F, No.18, Ln. 568, Niupu E. Rd., Xiangshan Dist.,
 Hsinchu 30091, Taiwan(R.O.C)
 TEL : +886-3-5308665
 FAX : +886-3-5308663

2. COMPOSITION, INFORMATION ON INGREDIENTS

Water: >50%(vv)	EC#231-791-2 CAS#7732-18-5
Glycerol: 10-25 % (v/v)	EC#200-289-5 CAS#56-81-5
Sodium Dodecyl Sulfate: 1-3 %	EC#205-788-1 CAS#151-21-3
Urea: 3.5-6.5 M	EC#200-315-5 CAS#57-13-6
Tris- Phosphate: 1-2 %	EC#230-260-2 CAS#6992-39-8

Other components:

Components not listed here are not dangerous or their concentrations do not exceed the limits specified in the EU directive 1999/45/EC.

Additional information: For the wording of the listed risk phrases refer to section 16.

3. HAZARDS IDENTIFICATION

The preparation is not classified as dangerous according to Directive 1999/45/EC or Directive 67/548/EC, and is not labelled with danger symbols or indications of danger.

Hazards deriving from physical-chemical properties:

The data available do not support any physical or chemical hazard.

Human health hazards:

Inappropriate handling may cause formation of aerosols. Inhalation of aerosols may cause allergic type reactions in sensitized individuals. Prolonged skin contact may cause minor irritation.

Caution: Avoid inhalation, ingestion, and skin contacts.

Environmental hazards: The data available do not support any environmental hazard.

4. FIRST AID MEASURES

Ingestion:

If swallowed, rinse mouth and throat thoroughly with water. Consult a doctor in event of any complaints.

Inhalation:

If inhaled, remove victim from area of exposure. If symptoms of irritation or respiratory allergy occur, consult a doctor.

Skin contact:

Generally, the product does not irritate the skin. If symptoms of irritation after prolonged exposure occur, consult a doctor.

Eye contact:

In case of contact, wash out opened eyes with plenty of water for several minutes. Consult a doctor in event of any complaints.

5. FIRE FIGHTING MEASURES

The preparation cannot cause fire.

Suitable extinguishing agents:

Water spray, Carbon Dioxide, dry chemical powder or appropriate foam.

Not suitable extinguishing agents: NONE

Special exposure hazards from decomposition products:

Not thoroughly investigated.

Protective equipment:

Wear self-contained breathing apparatus and protective clothing

6. ACCIDENTAL RELEASE MEASURES

Personal precautions:

Wear appropriate gloves (see Section 8), safety glasses to clean small releases/spillage. In case of vapor/aerosol formation, wear suitable respiratory protection means. Ventilate the area.

Environmental precautions: No special measures required.

Methods for cleaning up:

Absorb spillages with liquid-binding material. Ventilate and wash the area after the spillage is completely removed.

Additional information:

No dangerous substances are released.

7. HANDLING AND STORAGE

Handling: No further requirements.

Use appropriate personal protection measures (see Section 8).

Avoid any contact with the liquid. Avoid the formation of aerosols. Ventilate area.

Storage:

Store at - 20°C.

None other specific requirements relating to the preparation.

8. EXPOSURE CONTROLS, PERSONAL PROTECTION

Ingredients with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that need to be monitored at the workplace.

General protective and hygienic measures:

The usual precautionary measures are the same as for chemicals handling.

Engineering measures:

None specific requirements relating to the preparation.

Personal protection:

- Respiratory protection: No special measures required.

- Hand protection: Wear rubber gloves. None special requirements to rubber.

- Eye protection: No need to wear glasses (goggles) when ordinary operations

with the preparation. Safety glasses are recommended during certain operations when aerosol formation is possible.

- Skin protection: The usual precautionary measures accepted in research laboratories, e.g. protective clothes (gowns).

None specific requirements

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Liquid.

Color: Colorless.

Odor: Odorless.

pH-value at 25°C: 8.0.

Boiling point: 100°C.

Flash Point: Not applicable.

Explosive properties: The preparation does not possess explosive properties.

Oxidizing properties: Not applicable.

Vapor pressure: Not applicable.

Relative density at 20°C: ~1 g/cm³.

Solubility in/miscibility with water: Fully miscible.

Fat solubility: Not applicable.

Vapor density: Not applicable.

Evaporation rate: Not applicable.

Conductivity: Not applicable.

Viscosity: Not viscous.

Product name: AccuRuler RGB Prestained Protein Ladder**Catalog number : 02101-250**

Data updated: 01/13/2022

MaestroGen Inc.
3F, No.18, Ln. 568, Niupu E. Rd.,
Xiangshan Dist., Hsinchu 30091, Taiwan
+886-3-5308665

10. STABILITY AND REACTIVITY

Stability:

Stable under storage conditions (see Section 7).

Conditions to avoid:

No conditions which may cause dangerous reactions.

Possibility of hazardous reactions:

No dangerous reactions known.

Materials to avoid:

No materials which may cause dangerous reactions.

Hazardous decomposition products:

No dangerous decomposition products known.

Hazardous polymerization: Will not occur.

11. TOXICOLOGICAL INFORMATION

Routes of exposure:

Inhalation: Inhalation of aerosols may cause allergic type reactions in sensitized individuals.

Ingestion: Repeated ingestion or excessive dosage may cause irritation of the gastrointestinal tract.

Skin and eye contact: Generally, the product does not irritate the skin and eyes. Prolonged skin or eye contact may cause minor irritation.

Additional toxicological information:

When used and handled according to specifications, the product does not have any harmful effects.

The product is not listed in a priority list as foreseen under Council Regulation (EEC) No 793/93 on the evaluation and control of the risks of the existing substances.

12. Ecological Information

Ecotoxicity:

The preparation is considered not dangerous to the environment with respect to mobility, persistence and degradability, bioaccumulation, aquatic toxicity and other data relating to ecotoxicity.

Persistence and degradability:

The preparation is readily biodegradable.

Bioaccumulative potential:

It is considered that the preparation does not accumulate in biota, and does not pass through the food chain.

Other adverse effects: Not applicable.

13. DISPOSAL CONSIDERATIONS

The disposal of the product does not present a danger.

No special disposal methods required, except to be in accordance with current local environmental regulations.

14. TRANSPORT INFORMATION

Special precautions: None known.

Proper shipping name:

Not dangerous goods according to transport regulations.

National transport regulations:

No additional national transport regulations are known to the supplier.

IMDG: Non-hazardous for sea freight.

ADR: Non-hazardous for road transport.

RID: Non-hazardous for rail transport.

ICAO/IATA: Non-hazardous for air transport.

Other information:

The preparation may be stored in ice during the transportation.

15. REGULATORY INFORMATION

US:

Labelling according to US Regulations: Caution: Avoid contact and inhalation.

US Regulatory information: SARA LISTED: No.

EU:

Labelling according to EU Guidelines: The preparation is not dangerous within the meaning of Directives 67/548/EC, 1999/45/EC, and 88/379/ECC, and is not labelled with danger symbols or indications of danger.

Canada:

WHMIS classification: This product has been classified in accordance with the hazard criteria of the CPR, and the MSDS contains all the information required by the CPR.

DSL: No.

NDSL: No.

16. OTHER INFORMATION

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. MaestroGen shall not be held liable for any damage resulting from handling or from contact with the above product.